



**ALT-1200 - Alerta Commode & Transfer Chair** - comes with drop down armrests, swing away, detachable and fold-up footrests, arm pads and rear brakes, with wipeable PVC seat and back. Seat is 53cm deep x 46cm wide.

### KEY FEATURES AND BENEFITS

- > Stable, robust powder coated steel frame
- > Drop down padded armrests
- > Swing away detachable footrests
- > Full length arm pads
- > Push handle for optimal positioning
- > 4 solid steel 5" swivel castors
- > Fold-up footrests
- > Wipeable PVC seat and back
- > Max user weight: 110kg