



Oak



Walnut

The **Alerta Gas Lift Overbed Table** features an easy lift handle to adjust the table height with minimal effort. It is finished in oak or walnut and comes with brakes on the two rear castors. The table has a laminate top and comes with a wheelchair accessible base design.

TECHNICAL SPECIFICATION

- > Wheelchair access base
- > Braking castors
- > Max weight capacity: 15kg
- > Height adjustment range: 80-105cm
- > Table top dimensions: 80 x 40cm

PRODUCT CODES

- ALT-GL-O - Alerta Gas Lift Overbed Table, Oak**
- ALT-GL-W - Alerta Gas Lift Overbed Table, Walnut**

KEY FEATURES AND BENEFITS



Rear castors



Easy lift height adjustment



Won't topple over even with 15kg on top!