



accentu8  
Comfort made infinite

## haven

Retro styled design, beautifully tailored, haven is elegantly contoured, compact and practical in every way. The design features the accentu8 'ripple' cushion back as a modular, flexible and supportive solution to the traditional 'waterfall' back. haven evokes a sense of nostalgia with a classic Scandinavian touch for a timeless modern look and a modern feel.

- ❖ individual removable modular back cushion pads
- ❖ easy grip arms to assist getting in, up and out of chair
- ❖ soft-support for stability and comfort
- ❖ removable duvet seat pad for cleaning or updating
- ❖ our unique 'swift-fit' easy chair assembly and disassembly
- ❖ clean-lined zipped back cover with easy access to motor
- ❖ integrated pressure relief seat pad options
- ❖ standard castors for easy manoeuvrability
- ❖ accessories – lumbar cushions, 10 day battery backup, USB port
- ❖ FSC® available on request – FSC licence code FSC-C113594



haven pictured – duo rise-recline function with power, shown in healthcare fabrics: Camira Poplar HPP02 lavender and Aspen HPA02 lavender

## built for purpose

Our technical standards and our products have the status of class 1 medical devices as defined by EU Directive 93/42 and we are CE registered with the MHRA (Medicines and Healthcare Regulatory Authority.)

All our products comply with the Furniture and Furnishings (Fire) (Safety) Regulations for general domestic use and contract use to crib 5.

### frame

Engineered and built to meet high contract performance standards in healthcare situations and high domestic use. Constructed to BS EN 16139: 2013 using FSC® timber. Birch ply and birch hardwood used for shape and stress rails glued and stapled. Screws and bolts used for load-bearing. Frame parts are CNC milled for accuracy of cut parts, consistency of assembly and a robust stable frame.

### seat and back suspension

The seat suspension is constructed from steel zig zag springs. The back suspension is elasticated webbing. Both seat, back and arms are padded with CMHR foam.

### back cushion

Standard specification is CMHR foam core with 100% polyester Dacron wrap for loft and comfort

### seat cushion

Standard seat duvet pad specification is CMHR foam core encased in fabric cover. The duvet wraps over the top of the seat box CMHR foam designed for support and comfort.

### mechanisms

German engineered mechanism (TUV type approved) made from high quality steel with a powder coated finish for durability

## our guarantee

5  
years

This model carries a five year warranty against manufacturing defects covering mechanisms and frames and one year on all fabrics. Subject to the terms and conditions.

## applications

This model is designed to be suitable for use in:

- ❖ hospitals
- ❖ care homes
- ❖ nursing homes
- ❖ residential retirement homes
- ❖ medical centres
- ❖ hospices
- ❖ clinics
- ❖ domestic homes

## fabrics

This model is available in a range of fabrics from our 'cover and colour' ranges suitable for contract and domestic use

### lifestyle 1 & 2

suitable for both domestic and residential retirement interiors and complement some general healthcare environments.

### healthcare

Specifically designed for high performance demands within healthcare offering anti MRSA, antimicrobial, antibacterial, antifungal, water repellent, bleach cleanable and stain resistance functionality. Many fabrics are FAG49 1000 NHS Contract approved

## environmental

accentu8 is an ethical, socially responsible company, Our furniture is made from certified sustainable timbers under the auspices of the Forest Stewardship Council® (FSC®). FSC® mix 70% (BV-COC-060268)

## product testing

All fabrics conform to BS5852 parts 1 and 2 test procedures for fire resistant levels 0 & 1 crib and match; wood crib 5 as in standard fabric.

Produced to BS7176 low and medium hazard areas of usage for both domestic and residential care

# select the style, select the function

With accentu8 there is no longer a need to compromise your choice with uncoordinated chair styles, comfort and functions. Our full range of power and manual mechanisms and functional chair operations are available across most models. We have made selecting user weights and seat sizes simple, helping you to make the right choice to meet the need.



## chair functions

### tilt-in-space rise-recline

- ❖ Controlled support, fixed position recline
- ❖ Angle of the hip and back remain constant
- ❖ Semi recline position, high leg lift
- ❖ Rises and assists from seated position
- ❖ Available as manual (recline only)
- ❖ Available as power one motor 2-way

### solo rise-recline

- ❖ Simple easy recline and relax operation
- ❖ chair reclines after footrest rises
- ❖ Back follows into recline position
- ❖ Rises and assists from seated position
- ❖ Available as manual (recline only)
- ❖ Available as power one motor 3-way

### duo rise-recline

- ❖ More power for more function
- ❖ Full recline and multi-positioning
- ❖ Footrest and back operate independently
- ❖ Rises and assists from seated position
- ❖ Available as power two motor 3-way

### trio rise-recline

- ❖ More power for more function
- ❖ Full recline and multi-positioning
- ❖ Footrest and back operate independently
- ❖ Rises and assists from seated position
- ❖ Available as power three motor 3-way

### trio rise-recline bariatric

- ❖ Uncompromised operation for bariatric users
- ❖ User weights up to 35, 40 and 50 stone
- ❖ Full recline and multi-positioning
- ❖ Footrest and back operate independently
- ❖ Footrest motor resists excessive leg weight
- ❖ Rises and assists from seated position
- ❖ Available as power three motor 3-way

### extended dimensions



1800 mm in  
full recline

# simply choose

## seat dimensions for haven

chair functions	model size	power rise/recline	manual recline only	user weight	seat height from floo	seat width	seat depth
tilt-in-space rise-recline	compact	+	+	25 st	16"	18"	18"
	standard	+	+	25 st	18"	20"	20"
	large	+	+	25 st	20"	23"	20"
solo rise-recline	compact	+	+	19 st	16"	18"	18"
	standard	+	+	19 st	18"	20"	20"
	large	+	+	19 st	20"	23"	20"
duo rise-recline	compact	+		19 st	16"	18"	18"
	standard	+		19 st	18"	20"	20"
	large	+		19 st	20"	23"	20"
trio rise-recline	compact	+		25 st	16"	18"	18"
	standard	+		25 st	18"	20"	20"
	large	+		25 st	20"	23"	20"
trio rise-recline bariatric	compact	+		35/40/50	18"	27"	20 - 24"
	standard	+		35/40/50	18"	29"	20 - 24"
	large	+		35/40/50	18"	31"	20 - 24"
	grand	+		35/40/50	18"	33"	20 - 24"

All copyright, design rights and other intellectual property in accentu8 products and designs are and remain the property of Wilcare Wales Ltd.

Wilcare Wales Ltd is a member of ACID (Anti Copying In Design)

All details are correct at time of publishing. All dimensions are approximate. accentu8 reserves the right to modify designs without prior notice.

accentu8 is a registered trademark of Wilcare Wales Ltd all chair designs are copyright of Wilcare Wales Ltd.

accentu8  
vale business park  
cowbridge cf71 7pf  
wales, united kingdom  
t 01446 776111  
f 01446 772226

info@accentu8.net  
www.accentu8.net

### Accreditations



The mark of responsible forestry