

OLYMPIC 180 HOIST



Product Code:

HLA034

Description:

- The Olympic 180 Hoist benefits from the rigorous ergonomic design procedures and quality checks.
- It has been constructed from high quality powder coated steel that will provide years of reliable service.
- The hoist boasts an impressive 180kg lifting capacity, powered leg opening and market leading Linak® actuators and motors.

Safe Working Load:

180kg (28 Stone)

Lifting Range:

1755mm

Product Length:

1265mm

Product Height:

1335mm/2010mm (highest position/lowest position)

Leg Inside Width:

490mm/1075mm (closed/open)

Castor Size:

75mm/100mm (Front/Rear)

Carrier Bar Type

3 Hooked safety attachments for loop slings

Key Features:

- Cost effective solution without compromising on quality
- Robust powder coated steel frame
- Linak control system
- Electric leg opening as standard
- Red emergency lowering screw
- Easy to position with 4 large swivel castors and ergonomic handle
- Tool free assembly
- Single person installation
- Can lift a patient from the floor

Warranty:

Subject to the terms of Harvest Healthcare Guarantee and Warranties the Hoist is guaranteed for five years from date of purchase.

Supplier Contact Details:

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk; w: www.harvesthealthcare.co.uk
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60 1BW

Important Safety Notice

Harvest Healthcare Ltd supplies equipment and appliances to cover a range of medical conditions. Specific choice of equipment is the sole responsibility of the clinical practitioner or carer overseeing individual patient care. Risk assessments specific to user must be carried out prior to the use of all products. Products and data continually undergo development and may therefore deviate from the details given.

Version 1: Jan 2021