

# HARVEST SLIDE SHEETS



Our slide sheets come in standard sizes to ease user movement and help carers.

Made from a low friction nylon, these slide sheets make the moving of the service user on a bed, or similar, an easy operation as well as assisting horizontal sides.

## CLEANING & CARE

For general cleaning, wipe with a soft cloth dampened with a mild detergent and water solution. This product may be machine washed as follows:

- 40C medium machine wash
- Cool iron
- Do not dry clean
- Do not tumble dry

## PLEASE NOTE

**In order to maintain this product's performance, it is advised that our slide sheets may only be machine washed a maximum of 20 times.**



**USE NON-SOLVENT AND NATURAL SOAP ONLY.  
DO NOT USE HYPERCARBONATE, PHENOL-BASED CLEANING  
SOLUTIONS, ABRASIVE COMPOUNDS OR CLEANING PADS.**