

## beautifully made built to last

A contemporary styled day chair, the Oxford is specially designed for healthcare and nursing care environments. Featuring large, easy to manoeuvre wheels with comfortable and sturdy push handles. Built in 7 day battery pack for cable free operation of the smooth tilt-in-space and extending foot rest action.

- mobile hoist friendly
- 2 motor tilt in space action
- 7 day remote power unit
- 4 hour quick charge
- Visco elastic pressure relief foam to seat
- independent operation of the footboard and backrest
- footrest folds under seat
- easy to manoeuvre 4 inch diameter castors (braked to rear)
- push handle to rear
- wide choice of fabrics
- optional wings
- optional lateral supports
- FSC® available on request – FSC licence code FSC-C113594
- optional seat cut out for third party



oxford

### \*sizes and mechanisms

model	user weight	seat width	seat depth	seat height	back height	overall width	overall depth
standard	19 stone	19"/21"	18"	18"	19"	32"	35"

model	Motors						Manual	
	1- 3 way	1TIS	2 - 3 way	2 TIS	3	2 + 4th position	3 way	TIS
standard				■				

The Wilcare Group of Companies,  
Vale Business Park Cowbridge  
CF71 7PF Wales, United Kingdom

† 01446 772227  
f 01446 772226  
sales@wilcarewales.co.uk  
www.wilcarewales.co.uk

Partner



\* 5 year warranty against manufacturing defects covers mechanisms and frames and one year on fabrics.  
\* Due to natural differences in types of fabrics and fillings variations in sizing can occur +/- 1/2  
Subject to the terms and conditions.