

## Bath Seat



The bath seat provides additional height and maximum comfort to individuals that find it difficult getting into and out of the bath, sitting on the base of the bath or have reduced movement.

The bath seat is constructed from sturdy, lightweight PVC plastic, with four rubber feet located on the bottom section which attach firmly to the bath to prevent slippage whilst in use.

## Specification

Product code	Dimensions
UBAB0006-2	460mmL x 280mmW x 150mmH
UBAB0008-2	460mmL x 280mmW x 200mmH
UBAB00012-2	460mmL x 280mmW x 300mmH

**Maximum user weight** 190 kg | 30 Stone

For more information on this product please contact us on  
tel: 0333 321 8996 or email: [sales@ultimatehealthcare.co.uk](mailto:sales@ultimatehealthcare.co.uk)

[www.ultimatehealthcare.co.uk](http://www.ultimatehealthcare.co.uk)

