

Wall Mounted Shower Stool

Product description

The wall mounted folding shower seat with legs, back and arms is intended to be used in a shower to provide a stable support, where a user is unable to stand for long periods and can be folded up against the wall when it is not in use.

The shower seat features hinged padded armrests on either side to provide excellent and safe support for the user when sitting or standing. The shower seat is fitted with a large, wide and deep back rest to improve user comfort, giving more support and reducing stress to the lower back.



Benefits

- Provides stable support for those who are unable to stand for long periods in the shower.
- Can be folded up against wall when not in use.
- Padded armrests to provide excellent and safe support when sitting or standing.
- Large, wide and deep back rest to improve user comfort.

Technical specification

Part number:	UBAS0040A	Max user weight:	190kg (30st)
Material:	Aluminium + Stainless Steel	Seat height:	450 – 575mm
Seat width:	460mm	Seat depth:	370mm



For more information on this product please contact us on **0333 321 8996** or email sales@ultimatehealthcare.co.uk

www.ultimatehealthcare.co.uk