

# Lifting Pole

Product Code: UBEB7003



The Lifting Pole and Handle provide stable overhead support for patients who find it difficult to rise out of bed.

Complimented by a moulded plastic height adjustable handgrip the lifting pole is designed for raising a patient to a sitting position comfortably. It can be fitted on the left or the right hand side of the bed..

## Specifications

**Overall height:** 1840mm

**Handle width:** 220mm

**Max user weight:** 75KG

### Ultimate Healthcare

Calmore Industrial Estate, Nutwood Way, Totton, Hampshire, SO40 3WW

tel: 0333 321 8996, fax: 0845 051 0035, email: [sales@ultimatehealthcare.co.uk](mailto:sales@ultimatehealthcare.co.uk)

[www.ultimatehealthcare.co.uk](http://www.ultimatehealthcare.co.uk)