

Integral Heel Support

Viscoelastic Foam with Memory Effect



The Integral heel support consists of viscoelastic foam with a memory effect which provides pressure reduction for those who are considered to be at a high risk of pressure ulcer development on their heels. The heel support is suited to users with normal or slightly limited mobility. It is also ideal for elderly users or those subject to constant movement of the lower limbs.

Features	Benefits	
Viscoelastic foam with memory effect	Viscoelastic foam gives the heel support better resistance against subsidence. The memory foam allows precise moulding of the heel therefore provides better pressure distribution.	
The Form	The Achilles tendon is protected by side supports and the forefoot and heel section are hollowed out to direct pressure towards lower risk regions.	
The Cover	The breathable, waterproof and vapour permeable cover promotes comfort whilst aiding infection control. The heel support is covered with a soft, bi-stretch polyurethane cover and has a comfortable attachment system offering sustainable support and easy installation.	

Technical Specification UHC Product code SYST'AM® Product code **Cushion Dimensions USYCOHEEL** P901T1HW 55 x 25x 7.5cm / 21¹/₂ x 10 x 3" Warranty: 2 years Cleaning of the cover: 55Cleaning of foam: Do not immerse in water

UHC product code: USYC0HEEL SYST'AM® product code: P901T1HW

Ultimate Healthcare | Calmore Industrial Estate | Nutwood Way | Totton | Hampshire | SO40 3WW tel: 0333 321 8996 | fax: 0845 051 0035 | email: sales@ultimatehealthcare.co.uk

www.ultimatehealthcare.co.uk