

Toilet Frame (Width Adjustable)

Product description

The height and width adjustable toilet frame is designed to provide support when toileting and is suitable for users who require help in standing or sitting.

This item is designed to be increased or decreased in both height and width. The toilet frame is strong, easy to use and ensures user comfort and safety.

The freestanding steel frame fits around the toilet and features moulded plastic arms which are attached to provide a comfortable grip and easy cleaning.



Benefits

- Provides support when toileting and is suitable for users who require help in standing or sitting.
- Height and Width adjustable.
- The toilet frame is strong, easy to use and ensures user comfort and safety.

Technical specification

Part number:	UTOF0021-2	Max user weight:	190kg
Overall height:	645-790mm	Overall width:	580mm – 730mm
Overall depth:	460mm		



For more information on this product please contact us on **0333 321 8996** or email sales@ultimatehealthcare.co.uk

www.ultimatehealthcare.co.uk