



BED CHANGING TROLLEY



Features & Benefits

Our Bed Changing Trolley is made with MRSA resistant coated steel and is designed for efficient bed changing. This trolley is designed to hold a single laundry bag and comes with a lid to help prevent any cross contamination between the clean and soiled linen as per CFPP management and provision guidelines. The trolley comes with three shelved trays which are removable for easy cleaning.

*Laundry bags sold separately



Specifications

STYLE	DESCRIPTION	SIZE
LDK/BCT	Bed Changing Trolley	102 x 108 x 47cm