

Product Specification



PRODUCT NAME: **Medical Pillow 29x52x5 cm- Royal Blue**

PRODUCT CODE: **3012057**

Size:	Width 52 cm, Length 29 cm, Height 5 cm
Filling:	Memory Foam density 45kg/m3 and size 50x26x5 cm, Soap shape
Fabric:	Polyurethane coated nylon, 180 g/m2, Royal Blue
Filter:	1S filter 9,5x27 cm
Label:	Screen print- white print on blue fabric
Pillow Weight:	660 g
Cleaning:	Disinfect with max. 72% alcohol, max. 1000 ppm active chlorine. Wipe off with a water moistened cloth
Safety Standards:	Flammability in accordance with BS EN ISO 12952-1:2010 and BS EN ISO 12952-2:2010

* Tolerances pillow weight +/- 3%; pillow size +/- 1cm

Gabriel Scientific OÜ

Mauri tee 3, Papsaare küla, Pärnu linn 88317, Estonia

Tel: +372 58119997

www.sleepangel-medical.com

SleepAngel.