



HMT106LO - CONNOISSEUR OVERBED TABLE



Features:

- Adjustable-height hand-operated steel column
- Height adjustment: 76 – 115cm
- Automatic safe lock
- Number of castors: 4 50mm (with brakes)
- Material: Frame: steel-epoxy.
- Vinyl wrapped trays

Recommendations:

- Do not lean heavily on the table because it could tip over.
- The table should not be used to support more than 10 kg.

Regular maintenance:

- Wipe away any spills on the shelf as moisture can damage the wooden tray.
- To clean, use a cloth lightly moistened with a soap solution, followed by drying.
- Do not use abrasive cleaners or bleach; this could damage or discolor the product.

Warranty:

For a period of 12 months from delivery faulty equipment will be repaired free of charge after receipt and examination by Harvest Healthcare [the Seller]. Paint finish and in-use damage together with operational negligence are not covered by warranty. The limit of the Seller's liability for warranty work is the credit of defective items returned to the Seller for examination. Unauthorised repairs or modifications to the equipment may nullify the warranty.



Harvest Healthcare Ltd

Sheaf House Bradmarsh Way
Bradmarsh Business Park
Rotherham S60 1BW. UK
+44(0)1709 377 172