

# Medi-roller Standard



**Medi-roller Standard transfer board** allows to move patients laterally comfortably and safely with less effort between beds, stretchers, operating tables, X-ray tables, etc.

Board made with an improved padded core more comfortable and light, covered with High Frequency welded polyamide fabric. Outer sliding sleeve made of very low friction silicone coated polyamide fabric. Also available disposable sleeve.

## MEDI-ROLLER STANDARD

### FEATURES



X-ray translucent and MRI safe.




Ultra sliding sleeve cover



High frequency welded board



Lightweight PE foam core



2 models: Large (foldable) Short

### SPECIFICATIONS

PRODUCT FAMILY: Medi-roller

PRODUCT NAME: Medi-roller Standard

#### SIZE CHART

Ref.	Model	Measures (cm)	Weight (Kg)
TMR6200	Foldable	180 x 50	2,00
TMR6220	Short	90 x 50	1,00

INTENDED USE: Lateral transfer of patients between beds, stretchers, operating tables, X-ray tables, etc.

BASIC UDI-DI: 843603535TMRNE

#### COLOUR

- Blue board
- Orange cover

#### COMPOSITION

**Board:** Ext: 70% PVC - 30% Polyester  
Int: Polyethylene

**Cover:** Silicone coated Polyamide  
LATEX FREE

#### LAUNDRY INSTRUCTIONS

**Board:** Clean with damp cloth and mild soap or surface disinfectant detergent.

**Cover:** - Machine wash max. 60°C - Do not tumble dry  
- Do not bleach - Do not iron  
- Do not dry clean




#### DISINFECTION RECOMMENDATIONS

- For infection prevention and control use the disposable protective cover (ref. TMR6270).
- Use a surface disinfectant or disinfectant wipes if required. For cover disinfection machine wash at 60°C using a disinfectant detergent. Do not use softener, chlorine nor whitening concentrates.

#### WARNING

- Maximum weight: 170 Kg.
- Maximum distance between surfaces: 20cm

#### PACKAGING



- 1 - Medi-roller Standard board
- 1 - Medi-roller Standard cover