

Toilet Aids



ALT-BE001

ALERTA TOILET SEAT AID, ADJUSTABLE HEIGHT

This toilet seat aid is plastic coated and has a large seating area. It features cut-outs at the front and rear for easy cleaning access. The toileting seat simply clicks into position on the frame and rotates when required. 165kg weight capacity.

- > Outer Product Width: 53-54cm
- > Outer Product Depth: 49.5-54cm
- > Overall Height: 69-83cm
- > Seat Height: 48-62cm
- > Boxed in 2's



ALT-BE003

ALERTA PORTABLE TOILET SURROUND, ADJUSTABLE HEIGHT

Portable toilet seat surround with height extension legs designed to make using the toilet easier for users in nursing homes and home care.

- > Outer Product Width: 55-60cm
- > Outer Product Depth: 49.5cm
- > Overall Height: 68-80.7cm
- > Boxed individually